ReSound GN

ReSound Assist: Kurzanleitung

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
The screenshot displays the ReSound Smart Fit 1.0 software interface. At the top right, a red box highlights the 'Online Services Login' button, with a circled '1' next to it. The interface is divided into several sections: 'Patient' (showing 'Michael Drew', 'First Time User', and 'Data Logging'), 'Instruments' (with 'Right' and 'Left' channels selected), 'Media Player', 'Initial Programs' (with dropdown menus for Program 1, 2, 3, and 4), and 'Right Instrument' and 'Left Instrument' details (including ReSound LiNX 3D 9, L7961-DRW, LP Receiver, Serial Number, Cable, and Battery). An audiogram graph is visible in the top right corner, showing hearing levels across frequencies from 125 to 8K Hz. The bottom right corner features a 'Fit Patient' button.

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
Patient Consent Needed

The sharing of health information needs patient consent every 12 months. Please read the document below and click continue.

Introduction

The purpose of this document is to allow you to take full advantage of all the features and functionalities in your hearing instruments manufactured by GN Hearing. To realize all the benefits from these advanced hearing instruments, your hearing care professional and GN Hearing must be able to receive and handle information related to your hearing loss via the GN Online Services.

When using the GN Online Services, information related to your hearing loss will be transmitted to enable direct communication, remote fine tuning and other services that enable you to hear more, do more and be more.

It is required by law that the manufacturer of your hearing instruments, GN Hearing, and your hearing care professional receive a consent from you prior to processing your personal information. All your information will be handled carefully, securely, and in accordance with the law. If you wish to review your personal information, or should you later wish to withdraw your consent, this can be done at any time. Details on information collection and purposes for processing are included below.

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Patient Consent Needed

The sharing of health information needs patient consent. Please check the box below and continue saving.

1

Accept and agree to the items identified, described and contained in SECTION ONE below

SECTION ONE

I consent to my Hearing Care Professional's disclosure to Manufacturer and onward transfer of the personal information in accordance with the description in the previous sections, including, but not limited to, types of personal information (incl. hearing impairment and other health information), purposes, recipients, etc. until I withdraw my consent to the disclosure or it lapses under applicable local laws:

I release my Hearing Care Professional from his/her obligation to maintain medical confidentiality within the boundaries of this consent form.

I agree to the processing of my personal information, including health information, and to the processing of uniquely identifying information (if any) for the following purposes and on the terms described in the "Privacy Information" section:

1) Processing and handling my remote request for remote fine tuning of my hearing instrument

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Patient Consent Summary

Please print and sign the document below. This should be kept in the Hearing Care Professionals records.

First Name: Last Name: Patient Email: 3

Introduction

The purpose of this document is to allow you to take full advantage of all the features and functionalities in your hearing instruments manufactured by GN Hearing. To realize all the benefits from these advanced hearing instruments, your hearing care professional and GN Hearing must be able to receive and handle information related to your hearing loss via the GN Online Services.

When using the GN Online Services, information related to your hearing loss will be transmitted to enable direct communication, remote fine tuning and other services that enable you to hear more, do more and be more.

It is required by law that the manufacturer of your hearing instruments, GN Hearing, and your hearing care professional receive a consent from you prior to processing your personal information. All your information will be handled carefully, securely, and in accordance with the law. If you wish to review your personal information, or should you later wish to withdraw your consent, this can be done at any time. Details on information collection and purposes for processing are included below.

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Request assistance [X]

BEFORE WE BEGIN

To check if we can help you with this service, please tell us if you're experiencing any of the following:

Has your hearing changed? Yes No
 Yes No

Is your hearing aid broken or clogged with wax/debris?
 Yes No

Is there a problem with the fit of the hearing aids in your ears?
 Yes No

Next [➔]

Request assistance [X]

DEFINE YOUR PROBLEM

Please check the problem(s) you are experiencing. Add your own description on the last page.

Too loud

Too soft

Muffled, unclear

Extra noise, whistling, static

Too sharp, tinny, or shrill

The sound fluctuates or cuts in and out

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Request assistance [X]

WHEN IS IT OCCURRING?

Please check when the problem occurs. Add your own description on the last page.

All the time

In certain environments

With specific sounds

Loud sounds: clanging dishes or cutlery, doors slamming, flushing water

Soft sounds such as keyboard clicks, footsteps, rustling paper

Speech, such as soft-spoken children when there is no background noise

My own voice

Wind noise

Previous [➔] Next [➔]

Request assistance [X]

HOW SEVERE IS THE PROBLEM?

Slightly annoying Annoying Hearing aid unusable

WHICH SIDE IS THE PROBLEM ON?

Left Both Right

WHICH PROGRAM(S)?

All programs

AllAround Restaurant

Previous [➔] Next [➔]

Request assistance [X]

HOW SEVERE IS THE PROBLEM?

Slightly annoying Annoying Hearing aid unusable

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
The screenshot displays the 'Fitting' screen of a hearing aid app. The top navigation bar includes 'Patient', 'Instruments', 'Fitting', and 'Summary'. The 'Fitting' tab is active, showing two frequency response graphs. The left graph shows a green line representing customer settings, and the right graph shows a grey line representing default settings. Below the graphs are two frequency response tables for 'All-Around' and 'Streaming Accessories' modes. The 'All-Around' table has a vertical slider set to '2'. At the bottom, there are buttons for 'Calibrate DFS', 'Pair Accessories', 'On', 'Off', 'View Patient Settings', and 'Save'. A red box highlights the 'View Patient Settings' button, with a callout '4' pointing to it. Below the main screen, a separate control bar shows a slider for 'Noise Tracker II' with settings 'Off', 'Mild', 'Moder...', 'Consid...', and 'Strong Per Env...'. A red box highlights the 'Strong Per Env...' setting, with a callout '2' pointing to it.

	All	250	500	750	1K	1.5K	2K	3K	4K	8K
All	0	8	10	12	18	20	21	20	18	
HL	0	3	5	7	7	10	12	11	10	
HD	0	1	0	1	1	4	6	4	4	
CR	1	2.2	2.8	2.8	2.8	2.3	2.1	2.2	1.7	
MPO	87	102	108	108	107	112	114	118	108	

	All	250	500	750	1K	1.5K	2K	3K	4K	8K
All	0	8	10	12	18	20	21	20	18	
HL	0	3	5	7	7	10	12	11	10	
HD	0	1	0	1	1	4	6	4	4	
CR	1	2.2	2.8	2.8	2.8	2.3	2.1	2.2	1.7	
MPO	87	102	108	108	107	112	114	118	108	

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GN Hearing GmbH
GN ReSound AG